independent Photographic Supplies Pty Ltd

'The Best Brands in Imaging'

ABN: 36 105 600 351
P.O. Box 317, Thornleigh NSW 2120

Telephone: +61 2 9875 2244 / Fax: +61 2 9875 5448
e-mail: sales@iphoto.net.au www.iphoto.net.au



IPS DRYLAB HEALTH CHECK

CHAMPION CSS Inks (Colour Specialist Series) Pre-installation Health Check for ALL Noritsu Drylabs*

For Noritsu Drylab *D701/ D703/ D1005/ QSS Green Simplex/ QSS Green Duplex & QSS Greenlab II Series Printers

Customer:	
Customer address:	
City/ Town:	State Post Code
IPS Customer Account No.	
Noritsu Drylab Model:	Serial No:
Test Print the following Images:	
Print Image 1: Printer's Nozzle Check	Does Print have broken lines? Yes/ No
Print Image 2: IPSSkyBeachTest.TIFF Check	Does Print have banding? Yes/ No
Print Image 3: IPSDryLabTestChart.TIFF	
Print Image 4: SureLabColorCheckSheet IPS.TIFF Print	If there are no broken line on Nozzle Check Prints
Print Image 5: IPS SkintoneAndBlack Test.TIFF	or Banding on "Skybeach test" Print:
Print Image 6: IPS BW Wedding Test.JPG	The Customer can order Champion CSS Inks
Customer and IPS Technical Rep to Sign & Date back of all abovementioned prints	
IPS Technical Rep to leave customer with copy of this document & copy of signed & dated prints	
Customer directed to: https://www.iphoto.net.au/champion.aspx for latest warranty document	
Return above images with original copy of this form to: Attn: Mr. Peter Hofland: Independent Photographic Supplies Pty Ltd P.O. Box 317,	
Thornleigh NSW 2120 AUSTRALIA	NOTE: PLEASE DO NOT FOLD IMAGES
This iPhoto Health Check Test was conducted on (DATE) :	
Customer's Name:	IPS Representatives Name:
Signature:	Signature:
Technotape CDCC (Felix Schoeller D) TETENOL Strakes EPSON	



















